Manual

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Medintim

A photograph of a man's bare torso. He is holding a manual pump in his right hand and a clear cylindrical container in his left hand. The pump is white and blue, with a blue handle. A blue line points from the pump's head to a text box. The background is a solid orange color.

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Manual

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A close-up photograph of a man's bare torso. He is holding a blue manual erection ring with both hands, positioning it around his waist. A blue line with a dot at the end points from the text box to the ring. The background is a solid orange color.

Medintim

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Manual

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A person is shown from the side, applying a blue, cylindrical device to their upper arm. The device has a white cap and a blue body. The person's arm is extended, and they are holding the device with their other hand. The background is a solid blue color.

Medintim

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A close-up photograph of a man's torso and arm. He is holding a clear, cylindrical medical device with a blue cap. He is using a small spray bottle to apply lubricant to the opening of the device. The background is a warm, orange-brown color.

Medintim

**Aplicați
gelul lubrifiant.**

Manual

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A photograph of a man's torso and legs, showing his midsection and thighs. He is holding a manual vacuum pump device. The device consists of a clear cylindrical chamber with a blue ring at the top and a grey handle with a blue lever. A blue line points from the text box to the opening of the cylinder where the penis is inserted.

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Manual

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A photograph of a man's torso and legs, showing he is using a manual vacuum pump. The pump is a clear cylindrical tube with a blue ring at the base and a grey handle with a blue trigger. The man's hands are visible, one holding the pump and the other adjusting the blue ring. The background is a solid orange color.

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